


## Foam Surfaces


### QM 3 — Zoned Preventative Surface

- Dual density polyurethane foam core
- Slight heel slope design 6" – 5"
- Strategic positioned contour cut top surface for optimum comfort and pressure reduction
- Contour cut reliefs in bottom surface for flexing to conform to bed deck
- 4-way stretch polyester cover
- Foam meets the fire retardant standards of California TB117
- Mattress meets 16 CFR1632 Flammability of Mattresses Standard
- Weight capacity: 350 lb.
- 36"x80"x6", Custom sizes available




### QM 4 — Zoned Preventative Surface with Side Bolsters

- Dual density polyurethane foam core
- Slight heel slope design 6" – 5"
- Firm raised bolster side rail 6" – 8" with gap for ingress/egress
- Strategic positioned contour cut top surface for optimum comfort and pressure reduction
- Contour cut reliefs in bottom surface for flexing to conform to bed deck
- 4-way stretch polyester cover
- Foam meets the fire retardant standards of California TB117
- Mattress meets 16 CFR1632 Flammability of Mattresses Standard
- Weight capacity: 350 lb.
- 36"x80"x6", Custom sizes available



### QM 6 — Premium Preventative Surface

- Multi density polyurethane foam core
- Slight heel slope design 6" – 5"
- Firm flat side bolster for edge protection
- Visco-elastic memory foam heel section
- Strategic positioned contour cut top surface for optimum comfort and pressure reduction
- Contour cut reliefs in bottom surface for flexing to conform to bed deck
- 4-way stretch polyester cover
- Foam meets the fire retardant standards of California TB117
- Mattress meets 16 CFR1632 Flammability of Mattresses Standard
- Weight capacity: 350 lb.
- 36"x80"x6", Custom sizes available



[www.quarthealthcare.com](http://www.quarthealthcare.com)

Ontario, Quebec & Atlantic Provinces Ph. 416.624.1431 e-mail: [pdunphy@quarthealthcare.com](mailto:pdunphy@quarthealthcare.com)  
BC, Alberta, Saskatchewan, Manitoba Ph. 778.239.9233 e-mail: [phil.mundy@quarthealthcarewest.com](mailto:phil.mundy@quarthealthcarewest.com)