

DynaLAL Low Air Loss Surface



The DynaLAL mattress replacement system combines 3-1 alternation and true low air loss features for the prevention and treatment of pressure ulcers.

With multiple pressure reduction features, simplicity of functional control, quick-disconnect components design, the DynaLAL Mattress Replacement System effectively manages ulcers and pain

NEW FEATURES:
Auto-Memory Reset
Lock-Out Capability

True Low Air loss Technology

72 openings (standard) on the air cells, constant air distributes equally under a vapour-permeable cover, wicking moisture away from the wound and creating a dry healing environment.

Deep Cell 3 to 1 Alternation Therapy

18 (standard) independent 8 inch high cells inflate and deflate in a 3 to 1 cycle, creating maximum support for the individual using 2/3 of the air cells.

Comfort Level

8 comfort levels for different weight ranges.

Static Mode

Provides constant low-pressure and acts as an alternative to the alternation mode.

Auto Firm

Quickly inflates the mattress with a simple push of a button for rapid installation and nursing care requirements.

Auto Alarm

Visual and audio alarms in the event of power failure.

CPR quick-Release Valve

Deflates the torso section in less than 30 seconds in the event of medical emergencies.

Mattress

Weight Capacity	Up to 1000 lbs
Inflated Dimension	36" x 80" x 9.5" (standard) *other widths and lengths available
Weight	33 lbs (standard)
Air Cell Material	Polyurethane
Cover Material	Quilted fabric with polyurethane coating
Base Material	Woven Polyester fabric with PVC backing
Safety Standard	Meets open flame requirements (FCC 16CFR 1633)

Control Unit

Inflation Time	Under 1 minute
Functions	Alternation Cycle (5 - 95 min) 8 comfort levels; Auto firm; Static mode Lock-out Capability, Auto-memory Reset Low Air Loss, CPR Valve
Power requirement	100~240 volt, 50/60Hz, 5 Amp
Size	16.5"(L) x 6.7"(W) x 10.6"(H)
Weight	13.2 lbs
Max. Air-flow	1200 lpm
Safety Standard	UL/cUL

www.quarthealthcare.com

Ontario, Quebec & Atlantic Provinces Ph. 416.624.1431 e-mail: pdunphy@quarthealthcare.com
BC, Alberta, Saskatchewan, Manitoba Ph. 778.239.9233 e-mail: phil.mundy@quarthealthcarewest.com