

## Turning System Lateral Rotation Low Air Loss Surface



The Turning System Lateral Rotation mattress replacement system is a low air loss mattress replacement system with lateral rotation therapy. It combines the most advanced technology and clinically effective features for the prevention and treatment of pressure ulcers.

With multiple pressure reduction features, simplicity of functional control, quick-disconnect components design, the Turning System Mattress Replacement System effectively manages ulcers and pain

### True Low Air loss Technology

With 72 openings (standard) on the air cells, constant air flow distributes equally under a vapour-permeable cover, wicking moisture away from the wound and creating an optimum healing environment.

### Rotation Angle

Can be independently selected for 20-degrees or 40-degrees. Assists with repositioning of individual.

### Comfort Level

8 comfort levels for different weight ranges. Permanently inflated side bolsters for added safety.

### Static Mode

Provides constant low-pressure and acts as an alternative to the rotation mode.

### Auto Firm

Quickly inflates the mattress with a simple push of a button for rapid installation and nursing care requirements.

### Auto Alarm

Visual and audio alarms in the event of power failure.

### CPR quick-Release Valve

Deflates the torso section in less than 20 seconds in the event of medical emergencies.

### Mattress

Weight Capacity	Up to 1000 lbs
Inflated Dimension	36" x 80" x 9.5"
Weight	33 lbs (standard)
Air Cell Material	Nylon
Cover Material	Quilted fabric with polyurethane coating
Base Material	Woven Polyester fabric with PVC backing
Safety Standard	Meets open flame requirements (FCC 16CFR 1633)

### Control Unit

Inflation Time	Under 1 minute
Functions	8 comfort levels; Auto firm; Static mode Lateral Rotation; Low Air Loss CPR Valve
Power requirement	100~120 volt, 60Hz, 5 Amp
Size	16.5"(L) x 6.7"(W) x 10.6"(H)
Weight	13.2 lbs
Max. Air-flow	1200 lpm
Safety Standard	UL/cUL

[www.quarthealthcare.com](http://www.quarthealthcare.com)

Ontario, Quebec & Atlantic Provinces Ph. 416.624.1431 e-mail: [pdunphy@quarthealthcare.com](mailto:pdunphy@quarthealthcare.com)  
BC, Alberta, Saskatchewan, Manitoba Ph. 778.239.9233 e-mail: [phil.mundy@quarthealthcarewest.com](mailto:phil.mundy@quarthealthcarewest.com)